



**GDP SWISS**  
*Beauty*  
PRODUKTKATALOG

i n s p i r i n g   c u s t o m e r s   w o r l d w i d e



## Introducing your New Practice Assistant



Cryo**Pen**<sup>®</sup>



**D**oes this sound familiar ? A patient comes into your practice for treatment of a cold, but that's not her biggest health problem. At a glance you can see a spot with actinic keratoses on the temple. The nature of the skin defect may require treatment. But will you bring up the subject ?

CryoPen offers the solution. The knowledge that in your practice you have the option of an instrument that can handle this treatment extremely accurately, quickly, effectively and with no stress, will allow you to address this imperfection to the patient's attention. With the CryoPen this procedure fits right into the flow of work in your practice.

**T**he Cryopen was designed with an obsession for details and a touch of originality. It brings out the perfectionist in everyone : check out the selection of applicators allowing for the treatment from a pinpoint 1mm up to a width of 20mm and more. From now you'll only treat the abnormal tissue with no more collateral damage to the healthy skin.



**More Quality**

# Accurate, practical

Cryosurgery, of course, is nothing new. It is recognized as THE reference in dermatology clinic. But to any doctors who've struggled with disposable cans of freezing spray, cryoliquid spills, imprecise cotton buds and laborious deliveries of hazardous liquid nitrogen, the CryoPen is a different proposition all together.



● Ø 1-3mm   ○ Ø 2-6mm   ● Ø 4-8mm   ● Ø 7-20mm

CryoPen®



The CryoPen emits a fine jet of nitrous oxide under high pressure, which allows you to work with millimeter precision. Now, there will be no more collateral damage to healthy tissue. It's so incredibly accurate and practical. You just load the device, choose the relevant applicator, point and freeze...

# and very safe...

**C**ryoPen's incredible precision delicately treats even the smallest skin defect in seconds. With no damage to healthy tissue, if necessary, your patient can now tolerate longer treatment, allowing for deeper freezing and enhanced outcomes. Requirement for repeated application will reduce dramatically and there is no problem with patient compliance, allergic reaction, irritation or damage to the healthy tissue around the lesion.

The CryoPen is a brilliant device for removing skin lesions quickly and effectively. Penetration of the freezing is 5 seconds for 1 millimeter. Consequently procedures typically last 5-25 seconds, so patients can be in and out of your clinic in just 10 minutes.



**S**ign up to watch the treatment videos on CryoPen treatment on : Angioma, Seborrheic Keratoses, Actinic Keratoses, Verrucae, Plantar Warts, Condyloma, Pigmented Spots, Solar Lentigo and others: [info@ho-equipments.com](mailto:info@ho-equipments.com)



***Actinic keratoses***, width approx. 6mm, depth approx. 2mm. White applicator, application time 12 seconds.



***Verrucae on hand***, width approx. 3mm, depth approx. 2mm. White applicator, application time 12 seconds.



### CryoPen|b+

#### **Excellence in cryotherapy made affordable.**

With 23.5g cartridges, this instrument is well equipped for the treatment of verrucae in GP practice  
2 applicators : White Dot and White dot 60 mm.  
The longer 60mm applicator is especially designed for condyloma treatment.



### CryoPen|m

#### **Accuracy at your fingertips.**

The size of the very elegant CryoPen|m is coupled with very practical 8g cartridges.  
2 applicators: Blue and White Dot.



### CryoPen|o

#### **Sufficient power at hand.**

Equipped for the use of 16g cartridges. Nice to have when power is required.  
2 applicators: White and Green Dot.



### CryoPen|x

#### **Cryotherapy at its best.**

With the possibility to use both 8g and 16g cartridges CryoPen|x can treat all skin imperfections you will be confronted with in your clinic.  
4 applicators : Blue, White, Green and Yellow Dot.  
2 applicators 60mm : Blue and White Dot.

### CryoPen|x+

#### **Extremely versatile.**

This version of CryoPen|x has 2 additional longer applicators (60mm) especially for condyloma treatment.  
4 applicators : Blue, White, Green and Yellow Dot +  
2 applicators 60mm : Blue and White Dot.

● Ø 1-3mm   ○ Ø 2-6mm   ● Ø 4-8mm   ● Ø 7-20mm



Penetration of the freezing in tissue is 1mm per 5 seconds. Consequently, the treatment of a 2 mm deep lesion will use 10 seconds of constant gas flow.

**Cartridges 8g:** 100 seconds of constant gas flow.

Packaging in box with 24 cartridges.

**Cartridges 16g :** 200 seconds of constant gas flow

Packaging in box with 6 cartridges.

**Cartridges 23,5g :** 300 seconds of constant gas flow

Packaging in box with 6 cartridges.





CryoPen is your direct source for diagnostic and therapeutic Health Care solutions. We invest a great deal of time and effort into research and development to maintain our technological lead.

#### Address

1, rue des Journaliers  
B - 7822 Ghislenghien

#### Telephone

0032 68 26 86 00

#### Email

[info@ho-equipments.com](mailto:info@ho-equipments.com)

[www.ho-equipments.com](http://www.ho-equipments.com)

---

**CDP SWISS**  
*Beauty*  
PRODUKTKATALOG

*Hauptsitz*

**CDP SWISS AG**  
Weinbergstrasse 22  
CH-8001 Zürich

*Niederlassung*

**CDP SWISS AG**  
Moosstrasse 31  
CH-8907 Wettswil am Albis

Tel. +41 43 538 57 45  
Fax +41 43 538 57 73

[www.cdpswiss.com](http://www.cdpswiss.com)